


Dameca PA230 VER

Vertical wall panel in high quality and flexibility



Dameca PA 230 VER is the ultimate wall panel with many opportunities for customised solutions that meet the growing demands of modern hospital buildings.

To ensure its user-friendly, functional and forward-thinking design, Dameca PA 230 has been developed in close cooperation with experienced industrial designers.

SPECIFICATIONS

- Flexible length from 0.5 m
- Aluminium construction with white plastic joints
- Available in spray-paint aluminium e.g. white
- Integrated gas outlets
- Easily accessible sockets, e.g. for data transfer and electricity
- Integrated patient alarm
- Easy to mount and move
- Fast and easy adaptation
- Optimal hygiene due to easy cleaning
- Maintenance-free
- Robust and damage-resistant

COMBINE THE PA230 VER WITH DUOFLEX RAIL

The Duoflex rail gives numerous possibilities for customised solutions. The rail can be mounted directly on the wall and be used as fitting for mounting e.g. TV/PC/entertainment system, oxygen flow meter and vacuum units.

The Duoflex rail is available in length from 0.30 m.

